

LED Stomacher

World's first homogenization program, fully automatic pedal distance system

Characteristics:

- World's first homogenization program.
- Fully automatic pedal distance system
- Larger and wider pedal size
- Capable of homogenizing 1,250ml/400ml
- Pedal detachable
- Close and open the door easily



Main Technical Parameters:

Model	TP-1250S	TP-400S
Timer(Unit)	1 second to 2 hours (sec) or Unlimited	
Blending speed	Steps 1 to 10 : 4.5 Strokes/sec ~ 12 Strokes/sec	
Control & Setup	Digital Control with Touch Screen	
LED Lamp	4 colors (Green:Closed, Blue:Open, White:Operating, Red:Waring)	
Power	Door and pedal completely detachable / interior washable	
Washable	24V high performance, low noise motor	
Motor	Hard, Medium, and Soft Auto programs depending on the food sample solidity	
Program	Dual Volt 110V / 220V	
Protection	Application of precision pedal distance sensor	
	Motor overload protection function	
	Motor Speed maintenance system applied	
Material	SUS 304 stainless	
Volume	1,250ml	400ml
LCD touch panel	7inch	5inch
Pedal distance	0 ~ 50(mm)	0 ~ 30(mm)
Weight	25kg	18kg
Size	(W)377 × (D)532 × (H)417mm	(W)280 × (D)440 × (H)260mm